

Product Datasheet

TBPL1 antibody (FITC) (orb350738)

Description	TBPL1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human TATA box-binding protein-like protein 1 protein (1-185AA)
Target	TBPL1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62380
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.