

Product Datasheet

Tpbg antibody (orb350768)

Description	Tpbg antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Rat Trophoblast glycoprotein protein (32-361AA)
Target	Tpbg
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5PQV5
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.