



Product Datasheet

ART5 antibody (orb35090)


Description	ART5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Ecto-ADP-ribosyltransferase 5 protein (45-291AA)
Target	ART5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96L15
Dilution Range	IHC-P:1:100-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat(lane 1),...



Immunohistochemical staining of human pa...



Immunohistochemical staining of human br...