

Product Datasheet

HLA-E antibody (FITC) (orb350954)

Description	HLA-E antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human HLA class I histocompatibility antigen, alpha chain E protein (33-261AA)
Target	HLA-E
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P13747
Expiration Date	12 months from date of receipt.