




Product Datasheet

SYTL2 antibody (orb35096)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	SYTL2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	104930
Uniprot ID	Q9HCH5
NCBI	NP_996810.1, NP_996813.1, NP_115755.2, NP_001156423.1, NP_001156425.1, NP_116561.1, NP_996812.1, NP_001156424.1, NP_996811.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of SYTL2 antibody ...



Flow cytometric analysis of HeLa cells I...