

Product Datasheet

Bontoxilysin-C1 antibody (Biotin) (orb350987)

Description	Bontoxilysin-C1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Clostridium botulinum Botulinum neurotoxin type C1 protein (2-428AA)
Target	Botulinum neurotoxin type C1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P18640
Expiration Date	12 months from date of receipt.