

Product Datasheet

Pectinase antibody (orb351000)

Description	Pectinase antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Cryptomeria japonica Polygalacturonase protein (46-433AA)
Target	Polygalacturonase
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P43212
Expiration Date	12 months from date of receipt.