

Product Datasheet

wnt8 antibody (Biotin) (orb351007)

Description	wnt8 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Xenopus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Xenopus laevis Protein Wnt-8 protein (23-358AA)
Target	wnt8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P28026
Expiration Date	12 months from date of receipt.