

Product Datasheet

hup1 antibody (FITC) (orb351035)

Description	hup1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Streptomyces lividans DNA-binding protein HU 1 protein (1-93AA)
Target	hup1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0A3H6
Expiration Date	12 months from date of receipt.