

Product Datasheet

dsbA antibody (Biotin) (orb351088)

Description	dsbA antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Pseudomonas syringae pv. tomato Thiol:disulfide interchange protein DsbA protein (23-214AA)
Target	dsbA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	052376
Expiration Date	12 months from date of receipt.