

Product Datasheet

RBMY1F antibody (Biotin) (orb351112)

Description	RBMY1F antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human RNA-binding motif protein, Y chromosome, family 1 member F/J protein (1-496AA)
Target	RBMY1F
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15415
Expiration Date	12 months from date of receipt.