

Product Datasheet

ST3GAL2 antibody (FITC) (orb351119)

Description	ST3GAL2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human CMP-N-acetylneuraminate-beta-galactosamide-alpha-2,3-sialyltransferase 2 protein (28-210AA)
Target	ST3GAL2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16842
Expiration Date	12 months from date of receipt.