

Product Datasheet

MAD2L1 antibody (FITC) (orb351127)

Description	MAD2L1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Mitotic spindle assembly checkpoint protein MAD2A protein (2-205AA)
Target	MAD2L1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13257
Expiration Date	12 months from date of receipt.