



Product Datasheet

KLDC2 antibody (orb35114)

Description	Rabbit polyclonal antibody to KLDC2. The protein is widely expressed, with high levels in skeletal m...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	KLHDC2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	46099
Uniprot ID	Q9Y2U9
NCBI	NP_055130.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of in mouse NIH-3T...