



Product Datasheet

EBP antibody (orb351145)


Description	EBP antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human 3-beta-hydroxysteroid-Delta(8),Delta(7)-isomerase protein (2-230AA)
Target	EBP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15125
Dilution Range	WB:1:200-1000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of A549(Lane 1), 2...



Immunohistochemical staining of human br...



Immunofluorescence analysis of 293T cell...