

## Product Datasheet

### Nqo2 antibody (FITC) (orb351179)

|                        |                                                                                       |
|------------------------|---------------------------------------------------------------------------------------|
| <b>Description</b>     | Nqo2 antibody (FITC)                                                                  |
| <b>Species/Host</b>    | Rabbit                                                                                |
| <b>Reactivity</b>      | Rat                                                                                   |
| <b>Conjugation</b>     | FITC                                                                                  |
| <b>Immunogen</b>       | Recombinant Rat Ribosylidihydronicotinamide dehydrogenase [quinone] protein (1-231AA) |
| <b>Target</b>          | Nqo2                                                                                  |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                           |
| <b>Form/Appearance</b> | Liquid                                                                                |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                         |
| <b>Note</b>            | For research use only                                                                 |
| <b>Isotype</b>         | IgG                                                                                   |
| <b>Clonality</b>       | Polyclonal                                                                            |
| <b>Uniprot ID</b>      | <b>Q6AY80</b>                                                                         |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                       |