

Product Datasheet

Hmgcs1 antibody (Biotin) (orb351240)

Description	Hmgcs1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mouse Hydroxymethylglutaryl-CoA synthase, cytoplasmic protein (1-520AA)
Target	Hmgcs1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8JZK9
Expiration Date	12 months from date of receipt.