

Product Datasheet

AEN antibody (FITC) (orb351263)

Description	AEN antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Apoptosis-enhancing nuclease protein (117-325AA)
Target	AEN
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WTP8
Expiration Date	12 months from date of receipt.