


Product Datasheet

N6AMT2 antibody (orb351265)

Description	N6AMT2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human EEF1A lysine methyltransferase 1 protein (2-214AA)
Target	EEF1AKMT1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WVE0
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat(Lane 1),...