



Product Datasheet

OXNAD1 antibody (orb351297)

Description	OXNAD1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Oxidoreductase NAD-binding domain-containing protein 1 protein (18-312AA)
Target	OXNAD1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96HP4
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse large intestine...



Western blot analysis of HeLa whole cell...