

Product Datasheet

KIAA1191 antibody (FITC) (orb351347)

Description	KIAA1191 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Putative monooxygenase p33MONOX protein (32-305AA)
Target	KIAA1191
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96A73
Expiration Date	12 months from date of receipt.