



Product Datasheet

AVL9 antibody (orb35137)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	AVL9
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Concentration	0.5 mg/ml
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	71947
Uniprot ID	Q8NBF6
NCBI	NP_055875.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse Neuro-2a ...