


Product Datasheet

DCAF4 antibody (orb351393)

Description	DCAF4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human DDB1- and CUL4-associated factor 4 protein (1-495AA)
Target	DCAF4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WV16
Dilution Range	WB:1:200-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T (Lane 1), H...