

Product Datasheet

GNL3L antibody (Biotin) (orb351464)

Description	GNL3L antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Guanine nucleotide-binding protein-like 3-like protein (1-582AA)
Target	GNL3L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NVN8
Expiration Date	12 months from date of receipt.