

Product Datasheet

L3MBTL4 antibody (HRP) (orb35148)

Description	L3MBTL4 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Lethal(3)malignant brain tumor-like protein 4 protein (1-270AA)
Target	L3MBTL4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NA19
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.