

Product Datasheet

GDAP2 antibody (FITC) (orb351491)

Description	GDAP2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Ganglioside-induced differentiation-associated protein 2 protein (1-496AA)
Target	GDAP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NXN4
Expiration Date	12 months from date of receipt.