

Product Datasheet

UTP6 antibody (Biotin) (orb351496)

Description	UTP6 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human U3 small nucleolar RNA-associated protein 6 homolog protein (1-220AA)
Target	UTP6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NYH9
Expiration Date	12 months from date of receipt.