

Product Datasheet

DEGP12 antibody (FITC) (orb351515)

Description	DEGP12 antibody (FITC)
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	FITC
Immunogen	Recombinant Arabidopsis thaliana Putative protease Do-like 12, mitochondrial protein (25-499AA)
Target	DEGP12
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9LK70
Expiration Date	12 months from date of receipt.