

Product Datasheet

MIS12 antibody (FITC) (orb351535)

Description	MIS12 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Protein MIS12 homolog protein (1-305AA)
Target	MIS12
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H081
Expiration Date	12 months from date of receipt.