

Product Datasheet

MAP1LC3A antibody (HRP) (orb351554)

Description	MAP1LC3A antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Microtubule-associated proteins 1A/1B light chain 3A protein (1-120AA)
Target	MAP1LC3A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H492
Expiration Date	12 months from date of receipt.