

Product Datasheet

POLG2 antibody (FITC) (orb351567)

Description	POLG2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human DNA polymerase subunit gamma-2, mitochondrial protein (1-485AA)
Target	POLG2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UHN1
Expiration Date	12 months from date of receipt.