



Product Datasheet

ARHGAP25 antibody (orb35158)

Description	Rabbit polyclonal antibody to ARHGAP25. The protein encodes negative regulators of Rho GTPases, whic...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	ARHGAP25
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	73435
Uniprot ID	P42331
NCBI	NP_055697.1 , NP_001159748.1 , NP_001159749.1 , NP_001007232.2
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of in mouse liver

...