

Product Datasheet

SLC13A4 antibody (HRP) (orb351602)

Description	SLC13A4 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Solute carrier family 13 member 4 protein (134-273AA)
Target	SLC13A4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UKG4
Expiration Date	12 months from date of receipt.