



Product Datasheet

SRPRB antibody (orb351625)

Description	SRPRB antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Signal recognition particle receptor subunit beta protein (58-271AA)
Target	SRPRB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y5M8
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of HeLa cells...



Western blot analysis of HepG2 whole ce...