

Product Datasheet

CYS4 antibody (Biotin) (orb351692)

Description	CYS4 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Saccharomyces cerevisiae Cystathionine beta-synthase protein (1-345AA)
Target	CYS4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P32582
Expiration Date	12 months from date of receipt.