


Product Datasheet

NAGK antibody (orb351705)

Description	NAGK antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human N-acetyl-D-glucosamine kinase protein (2-344AA)
Target	NAGK
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UJ70
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 (Lane 1), ...