

Product Datasheet

Capsid protein C antibody (orb351713)

Description	Capsid protein C antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Hepatitis C virus genotype 1a Genome polyprotein protein (192-325AA)
Target	Genome polyprotein
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P26664
Expiration Date	12 months from date of receipt.