

Product Datasheet

AIFM1 antibody (FITC) (orb351731)

Description	AIFM1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Apoptosis-inducing factor 1, mitochondrial protein (103-612AA)
Target	Minor capsid protein VP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95831
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.