

Product Datasheet

IL1F10 antibody (Biotin) (orb351748)

Description	IL1F10 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Interleukin-1 family member 10 protein (1-152AA)
Target	IL1F10
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WWZ1
Expiration Date	12 months from date of receipt.