

Product Datasheet

ST6GALNAC3 antibody (Biotin) (orb35177)

Description	ST6GALNAC3 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 3 protein (40-305AA)
Target	ST6GALNAC3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NDV1
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.