


Product Datasheet

UBE2D3 antibody (orb351773)

Description	UBE2D3 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Ubiquitin-conjugating enzyme E2 D3 protein (1-147AA)
Target	UBE2D3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P61077
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...