

Product Datasheet

gluD antibody (FITC) (orb351815)

Description	gluD antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Clostridium difficile NAD-specific glutamate dehydrogenase protein (1-421AA)
Target	gluD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P27346
Expiration Date	12 months from date of receipt.