

Product Datasheet

ANXA2 antibody (FITC) (orb351835)

Description	ANXA2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Annexin A2 protein (2-339AA)
Target	Gamma-glutamyl hydrolase
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07355
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.