

Product Datasheet

MRPL1 antibody (FITC) (orb351843)

Description	MRPL1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human 39S ribosomal protein L1, mitochondrial protein (51-324AA)
Target	MRPL1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BYD6
Expiration Date	12 months from date of receipt.