

Product Datasheet

Sialic acid-binding lectin antibody (FITC) (orb351851)

Description	Sialic acid-binding lectin antibody (FITC)
Species/Host	Rabbit
Reactivity	Frog
Conjugation	FITC
Immunogen	Recombinant Rana japonica Sialic acid-binding lectin protein (1-111AA)
Target	Sialic acid-binding lectin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P18839
Expiration Date	12 months from date of receipt.