

Product Datasheet

lolA antibody (Biotin) (orb351872)

Description	lolA antibody (Biotin)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Escherichia coli O9:H4 Outer-membrane lipoprotein carrier protein (22-203AA)
Target	lolA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A7ZYJ5
Expiration Date	12 months from date of receipt.