

Product Datasheet

L2 antibody (orb351877)

Description	L2 antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human papillomavirus type 18 Minor capsid protein L2 protein (1-462AA)
Target	L2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P06793
Expiration Date	12 months from date of receipt.