

Product Datasheet

ampC antibody (Biotin) (orb351888)

Description	ampC antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Pseudomonas aeruginosa Beta-lactamase protein (27-397AA)
Target	ampC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P24735
Expiration Date	12 months from date of receipt.