

Product Datasheet

nikR antibody (orb351969)

Description	nikR antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Escherichia coli Nickel-responsive regulator protein (1-133AA)
Target	nikR
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0A6Z6
Expiration Date	12 months from date of receipt.