

Product Datasheet

CAPN2 antibody (FITC) (orb352019)

Description	CAPN2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Porcine
Conjugation	FITC
Immunogen	Recombinant Pig Calpain-2 catalytic subunit protein (1-324AA)
Target	CAPN2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P43367
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.